



## Product Data Sheet

---

Product Name:	Rabbit Polyclonal Antibody to Human Amnionless (Rb460)
Product Number:	AG-0016-200
Antigen:	Synthetic peptide from the C-terminus of human amnionless ( <a href="#">AMN</a> ) (accession number <a href="#">AAK28532</a> ) prepared as a MAP peptide (Multiple Antigenic Peptide).
Type:	Rabbit polyclonal
Form:	IgG, Protein G-affinity purified
Specificity:	Recognizes human amnionless
Application:	ELISA (WB and IHC not evaluated)
Species Reactivity:	Human
Titer:	1:50,000 (ELISA measured against the MAP peptide)
Buffer:	Tris buffered saline, pH 7.4
Long Term Storage	-20°C
Shipping:	Ship on Blue Ice
Reference:	This antiserum was generated to a synthetic peptide corresponding to residues 415-430 (ASEELPLPRRLSLVPK) in the human amnionless.
Background Information:	Amnionless ( <a href="#">AMN</a> ) is a type I transmembrane protein that binds to the amino-terminal third of cubilin and directs subcellular localization and endocytosis of cubilin ( <a href="#">Fyfe et al. 2003</a> ).

---

**WARNING:** THIS PRODUCT IS NOT INTENDED OR APPROVED FOR HUMAN, DIAGNOSTICS OR VETERINARY USE. USE OF THIS PRODUCT FOR HUMAN OR ANIMAL TESTING IS EXTREMELY HAZARDOUS AND MAY RESULT IN DISEASE, SEVERE INJURY, OR DEATH.

**MATERIAL SAFETY DATA:** This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Material Safety Data Sheet.

**WARRANTY AND LIMITATION OF REMEDY:** Array Genetics, LLC makes no warranty of any kind, expressed or implied, including, but not limited to, the warranties of fitness for a particular purpose and merchantability, which extends beyond the description of the chemicals on the face hereof, except that the material will meet our specifications at the time of delivery. Buyers exclusive remedy and the sole liability of Array Genetics, LLC hereunder shall be limited to refund of the purchase price of, or at the option of Array Genetics, LLC the replacement of, all material that does not meet our specifications. Array Genetics, LLC shall not be liable otherwise or for incidental or consequential damages, including, but not limited to, the costs of handling. Said refund or replacement is conditioned on Buyer giving written notice to Array Genetics, LLC within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within said thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

---